

# specifications

The medical gas outlet shall be an Amico Alert-1 series.

Outlets shall be manufactured with a 6 - 1/2" [165 mm] length type "K" 1/2" [12.7 mm] outside diameter (3/8" nominal) size copper inlet pipe stub which is silver brazed to the outlet body. The body shall be of 1 - 5/16" [33 mm] diameter, one piece brass construction. For positive pressure gas services, the outlet shall be equipped with a primary and secondary check valve and the secondary check valve shall be rated at a maximum of 200 psi [1,379 kPa] in the event the primary check valve is removed for maintenance. Outlet bodies shall be gas specific by indexing each gas service to a gas specific dual pin indexing arrangement on the respective identification module.

A large color coded front plate shall be used for ease of gas identification and aesthetic appeal.

With the back rough-in mounted, the outlet shall adjust up to 3/4" [19 mm] variation in mounting plate thickness.

The latch-valve assembly shall be Oxequip/Medstar type and only accept corresponding Oxequip/Medstar type gas specific adapters.

All M.R.I. outlets are manufactured from nonferrous materials.

All outlets shall be cleaned and degreased for medical gas service, factory assembled and tested.

Amico products comply with NFPA-99.



# features

- ▶ Flow rate exceeds the requirements of NFPA
- ▶ Accepts only Oxequip / Medstar gas specific adapters
- ▶ 100 % Hydrostatic testing of each outlet and cleaned for medical gas service
- ▶ Indexed to eliminate interchangeability of gas services
- ▶ Up to 3/4" [19 mm] mounting plate thickness adjustment
- ▶ 360° swivel inlet pipe for easy installation

**project**


---



---



---



<sup>1</sup> NOTE: For CSA/ISO, suction inlets are not supplied with a secondary check valve. A pressure plug is provided for testing purposes, rated at 150 psi [1,034 kPa].

## model numbers

Front Latch Valve:

**O-FASC-OX-L-XXX**

Console Assembly:

**O-OXCON-L-XXX**

The L defines the language:

- U = English (NFPA)\*  
Oxygen, Medical Air, and  
Medical Vacuum
- E = English (CSA/ISO)
- F = French (CSA/ISO)

The XXX defines the gas:

- OXY = Oxygen
- AIR = Medical Air
- VAC = MedVac<sup>1</sup>
- CO2 = Carbon Dioxide
- N2O = Nitrous Oxide
- WAG = Waste Anesthetic  
Gas Disposal (NFPA)
- AGS = Anesthetic Gas Scavenging  
System<sup>1</sup> (CSA)

<sup>1</sup>Please see note above

## M.R.I. model numbers

Example: English (NFPA) Console Assembly = O-OXCONM-U-OXY

Console Assembly:

**O-OXCONM-L-XXX**

Front Latch Valve:

**O-FAMC-OX-L-XXX**

represented by: